REUSABLE CABLES FOR LANTERN LARYNGEAL ELECTRODE**Intended use**

Reusable cables for Lantern Laryngeal Electrodes are intended to connect Neurosign® Lantern Laryngeal Electrodes to intraoperative neuromonitoring systems for monitoring of the recurrent laryngeal or vagus nerves during surgery. Reusable cables for Lantern Laryngeal Electrodes are accessories to Neurosign® Lantern Laryngeal Electrodes and intraoperative neuromonitoring systems, and as such the cables are intended to be used by a qualified healthcare professional. This product can be used on the general population at the discretion of a healthcare professional.

Use of product

- Before use, clean and disinfect cables according to the instructions below.
- The 4203-NM cable is designed to be used in combination with Neurosign® Lantern Laryngeal Electrodes and Neurosign® V4 Intraoperative Nerve Monitor. The colour coding and design of connectors of this cable allow it to be used with other intraoperative neuromonitoring systems which utilise the same colour coding.
- The 4203-00 cable is designed to be used in combination with Neurosign® Lantern Laryngeal Electrodes and Neurosign® 100 Nerve Monitoring System.
- Reusable cables have 5 connectors with 6 wires coloured red, blue and green on the 4203-NM cable and yellow, blue and green on the 4203-00. These should be connected to the appropriate coloured connectors on the pre-amplifier pod.

Cleaning

Reusable cables for Lantern Laryngeal Electrodes should be thoroughly cleaned before the first use and after each subsequent use. Immediately after use, wipe the cable using a soft cloth dipped in lukewarm water (at temperatures below 43°C (110°F)) with a mild detergent/enzymatic product to remove foreign material. Excessive pressure and continued flexing of the cable during cleaning or use can stretch the cable jacket and break the internal wires, damaging the cable. After wiping clean, rinse the cable with lukewarm water to prevent coagulation and thus assisting in the removal of protein substances. Carefully inspect the cable to ensure the complete removal of any contamination. Do not immerse or rinse cable connectors in any solution, as it may damage the product.

REUSABLE CABLES FOR LANTERN LARYNGEAL ELECTRODE

Note: Avoid contact with strong, aromatic, chlorinated, ketone, ether, or ester solvents. Prolonged immersion (more than 1 hour) in alcohol or mild solvents, detergent solutions or highly alkaline solutions may affect the cable and electrode jacket flexibility.

Disinfection

After cleaning, the reusable cable for Lantern Laryngeal Electrodes may be disinfected with the use of 70% ethanol. The cable should be wiped with a soft cloth dipped in a 70% ethanol for at least 1 minute, to ensure proper disinfection of the product.

Drying

After cleaning and disinfection, dry the cable with a clean, soft, low-linting cloth to remove any remaining water or disinfection solution.

Sterilization

Reusable cables for Lantern Laryngeal Electrodes are not designed to sustain any type of sterilization, as it may damage the product.

Inspection

Visually inspect the cable for deterioration or excessive wear prior to each use. If cable or connectors have any fracture, damaged wire insulation, or unnatural angle, discard the cable according to your hospital's procedure for medical waste. The integrity of the cable can be checked with a simple resistance meter.

Storage and Transport

Protect cables from damage and contamination during storage and transport. When not in use, cables should be loosely coiled and stored at room temperature. Do not wrap cables around the equipment casing to avoid damaging internal wires. When transporting the products outside their main area of use, pack them in a resealable plastic bag to prevent damage or contamination.

Note: All serious incidents associated with the use of reusable cables for Lantern Laryngeal Electrodes should be reported to Technomed Europe at quality@technomed.nl and to a competent authority of the country where user or patient is established.

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The Magstim Company Ltd.
ACRA Regulatory Services Pty Ltd