DMNE

DISPOSABLE MONOPOLAR NEEDLE ELECTRODES

Detachable or pre-wired:

- **Pencil point**
- **Trocar tip**
- **Sensory**

INSTRUCTIONS FOR USE

EN

DISPOSABLE MONOPOLAR NEEDLE ELECTRODES

Detachable or pre-wired: Pencil point, Trocar tip, Sensory

Intended use and indications for use

Disposable Monopolar Needle Electrodes are intended for use with recording and monitoring equipment for the recording of biopotential signals including electroencephalography (EEG), electromyography (EMG) and nerve potential signals, and are intended for stimulation and recording with stimulation / recording equipment for EMG and nerve potential signals.

These products are indicated to be used for the recording, stimulation and monitoring of EMG, EEG and nerve potential signals in patients undergoing clinical neurophysiological testing.

Intended users

The intended users are healthcare professionals specifically trained and certified in electrophysiology techniques.

Patient target groups

The patient target groups are any patients identified by these medical experts to benefit from the diagnostic procedures as described above.

Use of product

* Select a device of suitable length and diameter for the investigation being undertaken.

* Safe removal of protection sheath: twist and pull the protection sheath, while securely holding the hub.

* Do not remove electrodes from patient by pulling on wires.

WARNING

If used as stimulation electrode, prolonged stimulation outputs above 0.25W/cm² could result in thermal burns.

If used as recording electrode, the possible loss of electrocautery return path may result in high current density passing through the needle, possibly resulting in tissue heating and damage.

Applying excessive force to the device may cause it to bend or break. In case the device breaks, parts may remain in the patient and need to be surgically removed.

Cautions

* Federal law (USA) restricts this device to sale by or on the order of a physician.

* Do not use the device if the sterile packaging is opened or damaged or if the use-before-date has passed.

* Inspect the device after opening in a sterile field. Discard if damaged.

* For single-patient use only. Do not reuse or re-sterilize. Cleaning and re-sterilizing the device can affect its safety, performance and effectiveness and expose patients and users to unnecessary risks such as infection and transmissible disease.

* Sterility is guaranteed up to expiration date unless package is opened or damaged.

* The device has been sterilized using ethylene oxide gas. Some patients may be allergic to residuals of this gas.

* Do not use the device in strong electrical fields, as any induced electrical field may influence the stimulation output on the connected equipment (if applicable), may make the readings on the connected equipment unreliable and may result in localized tissue heating.

* Do not insert the device into patients all the way up to the hub.

* Do not straighten or re-insert the device if it bends before, during or after insertion. Bent devices should be discarded and replaced with a new item.

MRI safety information

Do not use in an MRI scanning environment. The device has not been evaluated for safety in an MR environment. It has not been tested for heating or unwanted movement in an MR environment. Its safety in the MR environment is unknown. Performing an MR exam on a person who has this medical device inserted or positioned on them may result in injury and / or device malfunction.

Compatibility with accompanying devices

This device is equipped with DIN 42802 connectors and is compatible with any electromyography equipment that is equipped with DIN 42802 connections. Variants with a detachable lead wire also require a separate reusable cable.

Disposal

Always discard used devices in a properly marked medical biohazard sharps container.

Reporting incidents

All serious incidents associated with the use of this device should be reported to Technomed Europe at quality@technomed.nl and to a competent authority of the country where the user is established.

Electronic instructions for use

Access the most recent electronic instructions for use through the link which is provided on the label of the device.

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Pty Ltd
The Magstim Company Ltd.