DSNE (Pt/Ir)

DISPOSABLE SUBDERMAL NEEDLE ELECTRODES (PLATINUM / IRIIDIUM)

STERILE EO

EN: INSTRUCTIONS FOR USE

BG:	ИНСТРУКЦИИ ЗА УПОТРЕБА
CS:	NÁVOD K POUŽITÍ
DA:	BRUGSANVISNING
DE:	GEBRAUCHSANWEISUNG
EL:	ΟΔΗΓΙΕΣ ΧΡΗΣΗΣ
ES:	INSTRUCCIONES DE USO
ET:	KASUTUSJUHEND
FI:	KÄYTTÖOHJEET
FR:	CONSIGNES D'UTILISATION
HR:	UPUTE ZA UPORABU
HU:	HASZNÁLATI ÚTMUTATÓ
IS:	NOTKUNARLEIÐBEININGAR
IT:	ISTRUZIONI PER L'USO
JA:	使用説明書
KO:	사용 설명서
LT:	NAUDOJIMO INSTRUKCIJOS
LV:	NORĀDĪJUMI PAR LIETOŠANU
MK:	УПАТСТВО ЗА УПОТРЕБА
NL:	GEBRUIKSAANWIJZING
NO:	BRUKSANVISNING
PL:	INSTRUKCJA STOSOWANIA
PT:	INSTRUÇÕES DE UTILIZAÇÃO
RO:	INSTRUCȚIUNI DE UTILIZARE
RU:	ИНСТРУКЦИЯ ПО ПРИМЕНЕНИЮ
SK:	POKyny NA POUŽITIE
SL:	NAVODILA ZA UPORABO
SR:	UPUTSTVO ZA UPOTREBU
SV:	BRUKSANVISNING
TR:	KULLANIM TALİMATLARI

DISPOSABLE SUBDERMAL NEEDLE ELECTRODES

(PLATINUM / IRIIDIUM)

Intended use and indications for use

Platinum / iridium Disposable Subdermal Needle Electrodes are intended to be inserted in the subdermal or muscle tissue to sense bio-electric (electromyographic [EMG] or electroencephalographic [EEG]) signals in the distal nerves, and to be connected to EMG / EEG recording equipment. Positioning and use takes place under supervision of a licensed physician.

These products are indicated for the recording and monitoring of EMG, EEG and nerve potential signals of peripheral, cranial and spinal nerve roots during clinical neurophysiology testing or during surgery combined with intraoperative neurophysiological monitoring.

Intended users

The intended users are healthcare professionals specifically trained and certified in electrophysiology techniques.

Patient target groups

The patient target groups are any patients identified by these medical experts to benefit from the diagnostic procedures as described above.

Use of product

- * Select a device of suitable length and diameter for the investigation being undertaken.
- * Safe removal of protection sheath: twist and pull the protection sheath, while securely holding the shrink tubing just below the connection point or by holding the hub.
- * Do not remove electrodes from patients by pulling on wire.

WARNING

This device is not intended for use as stimulating electrode or patient ground electrode.

If used as recording electrode, the possible loss of electrocautery return path may result in high current density passing through the needle, possibly resulting in tissue heating and damage.

Applying excessive force to the device may cause it to bend or break. In case the device breaks, parts may remain in the patient and need to be surgically removed.

Cautions

- * Federal law (USA) restricts this device to sale by or on the order of a physician.
- * Do not use the device if the sterile packaging is opened or damaged or if the use-before-date has passed.
- * Inspect the device after opening in a sterile field. Discard if damaged.
- * For single-patient use only. Do not reuse or re-sterilize. Cleaning and re-sterilizing the device can affect its safety, performance and effectiveness and expose patients and users to unnecessary risks such as infection and transmissible disease.
- * Sterility is guaranteed up to expiration date unless package is opened or damaged.
- * The device has been sterilized using ethylene oxide gas. Some patients may be allergic to residuals of this gas.
- * Do not use the device in strong electrical fields, as any induced electrical field may influence the stimulation output on the connected equipment (if applicable), may make the readings on the connected equipment unreliable and may result in localized tissue heating.
- * Do not insert the device into patients all the way up to the hub.
- * Do not straighten or re-insert the device if it bends before, during or after insertion. Bent devices should be discarded and replaced with a new item.

MRI safety information

Do not use in an MRI scanning environment. The device has not been evaluated for safety in an MR environment. It has not been tested for heating or unwanted movement in an MR environment. Its safety in the MR environment is unknown. Performing an MR exam on a person who has this medical device inserted or positioned on them may result in injury and / or device malfunction.

Compatibility with accompanying devices

This device is equipped with DIN 42802 connectors and is compatible with any intraoperative neuromonitoring equipment that is equipped with DIN 42802 connections.

Disposal

Always discard used devices in a properly marked medical biohazard sharps container.

Reporting incidents

All serious incidents associated with the use of this device should be reported to Technomed Europe at qualitytechnomed@welcony.com and to a competent authority of the country where the user is established.

Electronic instructions for use

Access the most recent electronic instructions for use through the link which is provided on the label of the device.

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Pty Ltd
The Magstim Company Ltd.