DCNE

DISPOSABLE (CONCENTRIC) NEEDLE ELECTRODES

Gold plated, Pt/Ir

STERILE EO

EN: INSTRUCTIONS FOR USE

BG:	ИНСТРУКЦИИ ЗА УПОТРЕБА
CS:	NÁVOD K POUŽITÍ
DA:	BRUGSANVISNING
DE:	GEBRAUCHSANWEISUNG
EL:	ΟΔΗΓΙΕΣ ΧΡΗΣΗΣ
ES:	INSTRUCCIONES DE USO
ET:	KASUTUSJUHEND
FI:	KÄYTTÖOHJEET
FR:	CONSIGNES D'UTILISATION
HR:	UPUTE ZA UPORABU
HU:	HASZNÁLATI ÚTMUTATÓ
IS:	NOTKUNARLEIÐBEININGAR
IT:	ISTRUZIONI PER L'USO
JA:	使用説明書
KO:	사용 설명서
LT:	NAUDOJIMO INSTRUKCIJOS
LV:	NORĀDĪJUMI PAR LIETOŠANU
MK:	УПАТСТВО ЗА УПОТРЕБА
NL:	GEbruiksaanwijzing
NO:	BRUKSANVISNING
PL:	INSTRUKCJA STOSOWANIA
PT:	INSTRUÇÕES DE UTILIZAÇÃO
RO:	INSTRUCȚIUNI DE UTILIZARE
RU:	ИНСТРУКЦИЯ ПО ПРИМЕНЕНИЮ
SK:	POKYN Y NA POUŽITIE
SL:	NAVODILA ZA UPORABO
SR:	UPUTSTVO ZA UPOTREBU
SV:	BRUKSANVISNING
TR:	KULLANIM TALİMATLARI
ZH-S:	使用说明

INSTRUCTIONS FOR USE

DISPOSABLE (CONCENTRIC) NEEDLE ELECTRODES Gold plated, Pt/Ir

EN

Intended use and indications for use

Disposable Concentric Needle Electrodes are intended to be inserted in the subdermal, muscle or nerve tissue to sense bio-electric, electromyography (EMG) or electroencephalography (EEG) signals distally and are intended to be proximally connected to EMG / EEG recording equipment.

These devices are indicated to be used for the recording of EMG and EEG signals in patients undergoing clinical neurophysiological testing.

Intended users

The intended users are healthcare professionals specifically trained and certified in electrophysiology techniques.

Patient target groups

The patient target groups are any patients identified by these medical experts to benefit from the diagnostic procedures as described above.

Use of product

- * Select a device of suitable length and diameter for the investigation being undertaken.
- * Safe removal of protection sheath: twist and pull the protection sheath, while securely holding the hub.
- * Do not remove electrodes from patients by pulling on wires.

WARNING

This device is not intended for use as stimulating electrode or patient ground electrode.

If used as recording electrode, the possible loss of electrocautery return path may result in high current density passing through the needle, possibly resulting in tissue heating and damage.

Applying excessive force to the device may cause it to bend or break. In case the device breaks, parts may remain in the patient and need to be surgically removed.

Cautions

- * Federal law (USA) restricts this device to sale by or on the order of a physician.
- * Do not use the device if the sterile packaging is opened or damaged or if the use-before-date has passed.
- * Inspect the device after opening in a sterile field. Discard if damaged.
- * For single-patient use only. Do not reuse or re-sterilize. Cleaning and re-sterilizing the device can affect its safety, performance and effectiveness and expose patients and users to unnecessary risks such as infection and transmissible disease.
- * Sterility is guaranteed up to expiration date unless package is opened or damaged.
- * The device has been sterilized using ethylene oxide gas. Some patients may be allergic to residuals of this gas.
- * Do not use the device in strong electrical fields, as any induced electrical field may influence the stimulation output on the connected equipment (if applicable), may make the readings on the connected equipment unreliable and may result in localized tissue heating.
- * Do not insert the device into patients all the way up to the hub.
- * Do not straighten or re-insert the device if it bends before, during or after insertion. Bent devices should be discarded and replaced with a new item.

MRI safety information

Do not use in an MRI scanning environment. The device has not been evaluated for safety in an MR environment. It has not been tested for heating or unwanted movement in an MR environment. Its safety in the MR environment is unknown. Performing an MR exam on a person who has this medical device inserted or positioned on them may result in injury and / or device malfunction.

Compatibility with accompanying devices

The device must be used with a reusable cable from the same brand having a 5-pole DIN connector.

Disposal

Always discard used devices in a properly marked medical biohazard sharps container.

Reporting incidents

All serious incidents associated with the use of this device should be reported to Technomed Europe at qualitytechnomed@welcony.com and to a competent authority of the country where the user is established.

Electronic instructions for use

Access the most recent electronic instructions for use through the link which is provided on the label of the device.

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Pty Ltd
The Magstim Company Ltd.