K-Wire Stimulation Lead, Sterile

STERILE EO

EN: INSTRUCTIONS FOR USE

BG:	ИНСТРУКЦИИ ЗА УПОТРЕБА
CS:	NÁVOD K POUŽITÍ
DA:	BRUGSANVISNING
DE:	GEBRAUCHSANWEISUNG
EL:	ΟΔΗΓΙΕΣ ΧΡΗΣΗΣ
ES:	INSTRUCCIONES DE USO
ET:	KASUTUSJUHEND
FI:	KÄYTTÖOHJEET
FR:	CONSIGNES D'UTILISATION
HR:	UPUTE ZA UPORABU
HU:	HASZNÁLATI ÚTMUTATÓ
IS:	NOTKUNARLEIÐBEININGAR
IT:	ISTRUZIONI PER L'USO
JA:	使用説明書
KO:	사용 설명서
LT:	NAUDOJIMO INSTRUKCIJOS
LV:	NORĀDĪJUMI PAR LIETOŠANU
MK:	УПАТСТВО ЗА УПОТРЕБА
NL:	GEBRUIKSAANWIJZING
NO:	BRUKSANVISNING
PL:	INSTRUKCJA STOSOWANIA
PT:	INSTRUÇÕES DE UTILIZAÇÃO
RO:	INSTRUCȚIUNI DE UTILIZARE
RU:	ИНСТРУКЦИЯ ПО ПРИМЕНЕНИЮ
SK:	POKYNÝ NA POUŽITIE
SL:	NAVODILA ZA UPORABO
SR:	UPUTSTVO ZA UPOTREBU
SV:	BRUKSANVISNING
TR:	KULLANIM TALİMATLARI

INSTRUCTIONS FOR USE

EN

K-WIRE STIMULATION LEAD, STERILE

Intended use and indications for use

Disposable K-Wire Stimulation Leads are intended to be used for the conduction of electrical signals, in combination with recording / stimulating equipment.

These products are applied for recording and stimulation during clinical neurophysiological testing and intraoperative neurophysiological monitoring to test, locate and identify peripheral, cranial and spinal nerve roots that are at risk of damage.

Intended users

The intended users are healthcare professionals specifically trained and certified in electrophysiology techniques.

Patient target groups

The patient target groups are any patients identified by these medical experts to benefit from the diagnostic procedures as described above.

Use of product

The device has one black wire. Connect the device to the blue or active connector on the stimulator pod. Use a separate needle electrode to act as the return.

WARNING

Careless use of any element in electrosurgical systems may cause serious burns. The use of the device with safely functioning and compatible electrosurgical generators and accessories is the user's responsibility. Read and understand all precautions and instructions before attempting to use any electrosurgical system.

When stimulating nerves, an appropriate voltage / current should be applied according to the stimulator manufacturer's operating instructions. If the applied current is too high, the device may get damaged.

Cautions

- * Federal law (USA) restricts this device to sale by or on the order of a physician.
- * Do not use the device if the sterile packaging is opened or damaged or if the use-before-date has passed.
- * Inspect the device after opening in a sterile field. Discard if damaged.
- * For single-patient use only. Do not reuse or re-sterilize. Cleaning and re-sterilizing the device can affect its safety, performance and effectiveness and expose patients and users to unnecessary risks such as infection and transmissible disease.
- * Sterility is guaranteed up to expiration date unless package is opened or damaged.
- * The device has been sterilized using ethylene oxide gas. Some patients may be allergic to residuals of this gas.
- * Do not use the device in strong electrical fields, as any induced electrical field may influence the stimulation output on the connected equipment (if applicable), may make the readings on the connected equipment unreliable and may result in localized tissue heating.
- * To avoid damage to the nerve monitor, the device must be kept separate from active electrosurgical devices when not in use. A separate holder is recommended.
- * Keep voltage / current settings at the lowest possible level to achieve the desired effects.
- * Device use and application is to be determined by the physician.

MRI safety information

Do not use in an MRI scanning environment. The device has not been evaluated for safety in an MR environment. It has not been tested for heating or unwanted movement in an MR environment. Its safety in the MR environment is unknown. Performing an MR exam on a person who has this medical device inserted or positioned on them may result in injury and / or device malfunction.

Storage, transport and handling

Keep away from sunlight.

Compatibility with accompanying devices

This device is equipped with DIN 42802 connectors and is compatible with any multi-channel nerve monitor or stimulation connector of standard electromyography system that is equipped with 1.5mm DIN 42802 connections. The other end of the device needs to be connected to a K-Wire, which is an implant with 2.4mm diameter.

To ensure proper functioning of the device, make sure the connector is fully inserted into the connection port on the equipment, and that it is correctly connected to the K-Wire.

Disposal

Always discard used devices in a properly marked medical biohazard container.

Reporting incidents

All serious incidents associated with the use of this device should be reported to Technomed Europe at qualitytechnomed@welcony.com and to a competent authority of the country where the user is established.

Electronic instructions for use

Access the most recent electronic instructions for use through the link which is provided on the label of the device.

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Pty Ltd
The Magstim Company Ltd.