

SUBDERMAL NEEDLE ELECTRODES



RELIABLE RESPONSE ACQUISITION FOR DYNAMIC ELECTROMYOGRAPHY RECORDING DURING SURGERY

- ✓ Packaged in sets to improve efficiency
- ✓ Suitable for use with multi-channel nerve monitors
- ✓ 14 mm subdermal needle
- ✓ 1.5 m twisted pair cables, which connect to the pre-amplifier pod for ease of use
- ✓ Cost-effective intraoperative nerve monitoring solutions

Neurosign® disposable subdermal needle electrodes are designed to acquire EMG signals from muscles.

SUBDERMAL NEEDLE ELECTRODES

2-CHANNEL ELECTRODE SET

PN 4469-00 (BOX OF 10 SETS)/5883-00 (SINGLE SET)

The 2-channel electrode set contains all required needle electrodes for 2-channel EMG monitoring, including a reference needle electrode (green) and a return needle electrode (white) for use with stimulating probes.

This electrode set includes:

- Blue electrodes, twisted pair (Ch1)
- Red electrodes, twisted pair (Ch2)
- Green electrode, single (Reference)
- White electrode, single (Return)

4-CHANNEL ELECTRODE SET

PN 4470-00 (BOX OF 10 SETS)/5882-00 (SINGLE SET)

The 4-channel electrode set contains all required needle electrodes for 4-channel EMG monitoring, including a reference needle electrode (green) and a return needle electrode (white) for use with stimulating probes.

This electrode set includes:

- Blue electrodes, twisted pair (Ch1)
- Red electrodes, twisted pair (Ch2)
- Violet electrodes, twisted pair (Ch3)
- Orange electrodes, twisted pair (Ch4)
- Green electrode, single (Reference)
- White electrode, single (Return)

8-CHANNEL ELECTRODE SET

PN 4472-00 (BOX OF 5 SETS)/5881-00 (SINGLE SET)

The 8-channel electrode set contains all required needle electrodes for 8-channel EMG monitoring, including a reference needle electrode (green) and a return needle electrode (white) for use with stimulating probes.

This electrode set includes:

- Blue electrodes, twisted pair (Ch1)
- Red electrodes, twisted pair (Ch2)
- Violet electrodes, twisted pair (Ch3)
- Orange electrodes, twisted pair (Ch4)
- Grey electrodes, twisted pair (Ch5)
- Yellow electrodes, twisted pair (Ch6)
- Brown electrodes, twisted pair (Ch7)
- Beige electrodes, twisted pair (Ch8)
- Green electrode, single (Reference)
- White electrode, single (Return)



Neurosign® is a registered trademark of The Magstim Company Limited. All specifications are subject to change. The information contained within this document is accurate at the time of publication.

© 2024 Technomed
N-2024-003.01_EN

Legal Manufacturer:
Technomed Europe
Amerikalaan 71
6199 AE Maastricht-Airport
The Netherlands

Tel: +31 43 408 68 68
infotechnomed@welcony.com

technomedmedical.com

Your NEUROSIGN® Business Partner: