

Disposable Monopolar EMG Needle Electrode – Pencil Point



Technomed understands the importance of patient comfort and EMG signal quality.

Technomed's Disposable Monopolar Needle Electrode – Pencil Point and cable solutions are designed for optimal balance between needle sharpness and ease of use.


The hub and cable provide a firm grip and optimal precision in needle positioning.

User Benefits

- Clear and reliable signal with improved recording area.
- Sharp needle tip for ease of penetration and repositioning.
- Approved for stimulation and recording.
- Precise handling with improved ergonomic hub and cable design.
- Optional reusable connection cable.

Available variants

Product Description	Needle Length	Diameter	Hub Color	Part Number
Pre-wired	25 mm	0.35 mm (28 Gauge)	Red	TE/AP2535-335
Pre-wired	25 mm	0.45 mm (26 Gauge)	Yellow	TE/AP2545-335
Pre-wired	37 mm	0.35 mm (28 Gauge)	Orange	TE/AP3735-335
Pre-wired	37 mm	0.45 mm (26 Gauge)	Green	TE/AP3745-335
Pre-wired	50 mm	0.45 mm (26 Gauge)	Blue	TE/AP5045-335
Pre-wired	75 mm	0.45 mm (26 Gauge)	Violet	TE/AP7545-335
Detachable	25 mm	0.35 mm (28 Gauge)	Red	TE/AD2535-008
Detachable	25 mm	0.45 mm (26 Gauge)	Yellow	TE/AD2545-008
Detachable	37 mm	0.35 mm (28 Gauge)	Orange	TE/AD3735-008
Detachable	37 mm	0.45 mm (26 Gauge)	Green	TE/AD3745-008
Detachable	50 mm	0.45 mm (26 Gauge)	Blue	TE/AD5045-008
Detachable	75 mm	0.45 mm (26 Gauge)	Violet	TE/AD7545-008




Packing: 25 individually packed electrodes per box
Pre-wired: 0.75 m

Reusable Cable	Wire Length	Part Number
Shielded	1.0 m	TE/D54-451
Unshielded	0.6 m	TE/D34-231

Specifications

Not made with natural rubber latex	✓
REACH & RoHS (EU) compliant	✓
EtO sterilized	✓
Lead wire core	Tin Plated Copper
Lead wire insulation	PVC
Connector	DIN 42802-1
Shelf life for unopened pouches	36 months



Your Technomed Business Partner:



Legal manufacturer:
Technomed Europe
 Amerikalaan 71

6199 AE | Maastricht-Airport | The Netherlands
 T +31 43 408 68 68 | info@technomed.nl

Technomed USA

1626 Terrace Drive

Roseville | MN 55113 | United States

T +1 651 288 0890 | info@technomed-usa.com

www.technomed.nl