

Disposable Concentric EMG Needle Electrode



Technomed understands the importance of a crystal clear signal, ease of use and patient comfort in EMG.

Technomed's Disposable Concentric Needle Electrode and cable are designed for optimal balance between needle sharpness and ease of use.


The hub and cable provide a firm grip and optimal precision in needle positioning.

User Benefits

- Clear signal with optimized recording area and shielded cable.
- Sharp tip for patient penetration and repositioning.
- Approved for recording.
- Precise handling with improved ergonomic hub and cable design.
- Color coded for ease of use.
- Easy positioning with recording area indicator mark.
- Reusable connection cable.

Available variants

Needle Length	Diameter	Recording Area	Hub Color	Part Number
25 mm	0.30 mm (30 Gauge)	0.031 mm ²	Red	TE/B50600-001
25 mm	0.40 mm (27 Gauge)	0.031 mm ²	Yellow	TE/B50600-002
30 mm	0.30 mm (30 Gauge)	0.031 mm ²	Pink	TE/B50600-006
37 mm	0.45 mm (26 Gauge)	0.085 mm ²	Green	TE/B50600-003
50 mm	0.45 mm (26 Gauge)	0.085 mm ²	Blue	TE/B50600-004
75 mm	0.60 mm (23 Gauge)	0.085 mm ²	Violet	TE/B50600-005



Packing: 25 individually packed electrodes per box

Specifications

Not made with natural rubber latex	✓
REACH & RoHS (EU) compliant	✓
EtO sterilized	✓
Shelf life for unopened pouches	36 months



Reusable Cable	Wire Length	Part Number
5-pole DIN connector	1.25 m	TE/D64-551
5-pole DIN connector	2.00 m	TE/D64-851



Your Technomed Business Partner:

USA
Rx ONLY



Legal manufacturer:

Technomed Europe

Amerikalaan 71

6199 AE | Maastricht-Airport | The Netherlands

T +31 43 408 68 68 | info@technomed.nl

Technomed USA

1626 Terrace Drive

Roseville | MN 55113 | United States

T +1 651 288 0890 | info@technomed-usa.com

www.technomed.nl