


# Disposable Hypodermic EMG Needle Electrode



**Technomed understands the importance of a crystal clear signal, ease of use, repositionability, and patient comfort.**

Technomed's Disposable Hypodermic Needle Electrode is a high quality product developed for EMG-guided drug delivery of Botulinum Toxin Type A.

Ergonomically designed hub and cable enables a firm grip and optimal precision in needle positioning.


#### **User Benefits**

- Sharp tip and special coating for skin penetration and repositioning.
- Leak-free and minimum drug residual volume.
- Optimized to fit all slip-tip and Luer Lock syringe connections.
- Color coded for ease of use.
- Single patient use to avoid cross-contamination.

## Available variants

Needle Length	Diameter	Lead Wire	Hub color	Part Number
25 mm	0.30 mm (30 Gauge)	0.75 m	Red	TE/L2530-335
25 mm	0.40 mm (27 Gauge)	0.75 m	Yellow	TE/L2540-335
37 mm	0.40 mm (27 Gauge)	0.75 m	Orange	TE/L3740-335
37 mm	0.45 mm (26 Gauge)	0.75 m	Green	TE/L3745-335
50 mm	0.50 mm (25 Gauge)	0.75 m	Blue	TE/L5050-335
75 mm	0.70 mm (22 Gauge)	0.75 m	Violet	TE/L7570-335

**Packing:** 10 individually packed electrodes per box



## Specifications

Not made with natural rubber latex	✓
REACH & RoHS (EU) compliant	✓
EtO sterilized	✓
Lead wire core	Tin Plated Copper
Lead wire insulation	PVC
Connector	DIN 42802-1
Shelf life for unopened pouches	36 months

Your Technomed Business Partner:



[www.technomed.nl](http://www.technomed.nl)

Legal manufacturer:

Technomed Europe

Amerikalaan 71

6199 AE | Maastricht-Airport | The Netherlands

T +31 43 408 68 68 | [info@technomed.nl](mailto:info@technomed.nl)

Technomed USA

1626 Terrace Drive

Roseville | MN 55113 | United States

T +1 651 288 0890 | [info@technomed-usa.com](mailto:info@technomed-usa.com)