

Reusable EEG Cup Electrode Gold and Silver



Technomed understands the importance of a crystal clear signal, ease of use and electrode durability.

Technomed's Reusable EEG Cup Electrode Gold and Silver are designed to provide reliable signals and longtime ease of use.

These electrodes are designed for easy adjustment to patient contours to avoid electrode lift.

The connector is optimized for easy headbox insertion and extraction.


User Benefits

- Clear and reliable signal.
- Short cup-to-cable connection to reduce lift risk.
- Optimal patient comfort with low profile disk, soft and flexible anti-tangle lead wires.
- Optimized connector design for easy headbox insertion and extraction.
- Long electrode life with watertight design and Kevlar reinforced wires.
- Available in two different cup sizes.

Available variants

Materials	Cup Size	Cup Height	Wire Length	Part Number
Silver Plated (Ag)	10 mm	3.1 mm	1.0 m	TE/C12-434
Silver Plated (Ag)	10 mm	3.1 mm	1.5 m	TE/C12-634
Silver Plated (Ag)	10 mm	3.1 mm	2.0 m	TE/C12-834
Silver Plated (Ag)	10 mm	3.1 mm	2.5 m	TE/C12-934
Silver Plated (Ag)	6 mm	1.7 mm	1.0 m	TE/C14-434
Silver Plated (Ag)	6 mm	1.7 mm	1.5 m	TE/C14-634
Silver Plated (Ag)	6 mm	1.7 mm	2.0 m	TE/C14-834
Gold Plated (Au)	10 mm	3.1 mm	1.0 m	TE/C32-434
Gold Plated (Au)	10 mm	3.1 mm	1.5 m	TE/C32-634
Gold Plated (Au)	10 mm	3.1 mm	2.0 m	TE/C32-834
Gold Plated (Au)	10 mm	3.1 mm	2.5 m	TE/C32-934
Gold Plated (Au)	6 mm	1.7 mm	1.0 m	TE/C34-434
Gold Plated (Au)	6 mm	1.7 mm	1.5 m	TE/C34-634
Gold Plated (Au)	6 mm	1.7 mm	2.0 m	TE/C34-834

Packing: sets of 12, 2 x 6 colors (red, green, white, yellow, blue, brown)



Specifications

Cup	Ag or Au plating
Not made with natural rubber latex	✓
REACH & RoHS (EU) compliant	✓
Lead wire design	Tinsel wire
Lead wire core	Kevlar and Silver Plated Copper
Lead wire insulation	TPU
Spring loaded connector	Gold Plated Brass / PVC
Connector	DIN 42802-1

Your Technomed Business Partner:

USA
Rx ONLY



Cleaning and disinfection: see IFU for details



www.technomed.nl

Legal manufacturer:

Technomed Europe

Amerikalaan 71

6199 AE | Maastricht-Airport | The Netherlands

T +31 43 408 68 68 | info@technomed.nl

Technomed USA

1626 Terrace Drive

Roseville | MN 55113 | United States

T +1 651 288 0890 | info@technomed-usa.com