


# Disposable EEG Cup Electrode Silver/Silver Chloride



## Technomed understands the importance of clear signals, efficiency and ease of use.

Technomed's Disposable Ag/AgCl EEG Cup Electrode is designed to optimize EEG procedures to enable focusing on the patient and the EEG recording.

The electrode is designed with flexible wires and optimized connector for easy headbox insertion and extraction.

## User Benefits

- No cleaning and preparation of EEG electrodes needed.
- Clear and reliable signal due to optimized layer of frequency stable Ag/AgCl .
- Short cup-to-cable connection to reduce lift risk.
- Optimized connector for easy headbox insertion and extraction.
- Optimized patient comfort with low profile disk, soft and flexible anti-tangle lead wires.
- Reduced risk of cross-contamination.


## Available variants

Materials	Cup Size	Cup Height	Wire Length	Packing	Part Number
Silver/Silver Chloride (Ag/AgCl)	10 mm	2.6 mm	1.0 m	Sets of 25, 5 x 5 colors	TE/C52-438
Silver/Silver Chloride (Ag/AgCl)	10 mm	2.6 mm	1.5 m	Sets of 25, 5 x 5 colors	TE/C52-638
Silver/Silver Chloride (Ag/AgCl)	10 mm	2.6 mm	2.0 m	Sets of 25, 5 x 5 colors	TE/C52-838
Silver/Silver Chloride (Ag/AgCl)	10 mm	2.6 mm	2.5 m	Sets of 25, 5 x 5 colors	TE/C52-938
Silver/Silver Chloride (Ag/AgCl)	10 mm	2.6 mm	1.0m	Sets of 10, 1 x 10 colors	TE/C52-432
Silver/Silver Chloride (Ag/AgCl)	10 mm	2.6 mm	1.5 m	Sets of 10, 1 x 10 colors	TE/C52-632
Silver/Silver Chloride (Ag/AgCl)	10 mm	2.6 mm	2.0 m	Sets of 10, 1 x 10 colors	TE/C52-832
Silver/Silver Chloride (Ag/AgCl)	10 mm	2.6 mm	2.5 m	Sets of 10, 1 x 10 colors	TE/C52-932
Silver/Silver Chloride (Ag/AgCl)	10 mm	2.6 mm	1.0 m	Sets or 27	TE/C52-4327
Silver/Silver Chloride (Ag/AgCl)	10 mm	2.6 mm	1.5m	Sets of 27	TE/C52-6327

**5 colors:** red, green, white, yellow, blue

**10 colors:** red, green, white, yellow, blue, grey, purple, orange, brown, black

**27 pack colors:** 6x red, 6x green, 6x white, 6x yellow, 6x blue



## Specifications

Cup	Ag/AgCl layer
Not made with natural rubber latex	√
REACH & RoHS (EU) compliant	√
Lead wire core	Tin Plated Copper
Lead wire insulation	PVC
Connector	DIN 42802-1

Your Technomed Business Partner:



USA
Rx ONLY



Legal manufacturer:

Technomed Europe

Amerikalaan 71

6199 AE | Maastricht-Airport | The Netherlands

T +31 43 408 68 68 | info@technomed.nl

Technomed USA

1626 Terrace Dr.,

Roseville | MN 55113 | United States

T +1 651 288 0890 | info@technomed-usa.com

[www.technomed.nl](http://www.technomed.nl)