

# Disposable Monopolar EMG Needle Electrode – Trocar Point



**Technomed understands the importance of a crystal clear signal, ease of use, and good grounding.**

**Technomed's Disposable Monopolar Needle Electrode – Trocar Point and cable solutions are designed for optimal sharpness and ease of use.**



**The hub and cable provide a firm grip and optimal precision in needle positioning.**

## **User Benefits**

- Clear and reliable signal with improved recording area.
- Sharp needle tip for optimized patient comfort, easy penetration and smooth repositioning.
- Approved for stimulation and recording.
- Precise handling with improved ergonomic hub and cable design.
- Optional reusable connection cable.

## Available variants

Product Description	Needle Length	Diameter	Hub Color	Part Number
Pre-wired	25 mm	0.35 mm (28 Gauge)	Red	TE/AP2535-335/01
Pre-wired	25 mm	0.45 mm (26 Gauge)	Yellow	TE/AP2545-335/01
Pre-wired	37 mm	0.35 mm (28 Gauge)	Orange	TE/AP3735-335/01
Pre-wired	37 mm	0.45 mm (26 Gauge)	Green	TE/AP3745-335/01
Pre-wired	50 mm	0.45 mm (26 Gauge)	Blue	TE/AP5045-335/01
Pre-wired	75 mm	0.45 mm (26 Gauge)	Violet	TE/AP7545-335/01
Detachable	25 mm	0.35 mm (28 Gauge)	Red	TE/AD2535-008/01
Detachable	25 mm	0.45 mm (26 Gauge)	Yellow	TE/AD2545-008/01
Detachable	37 mm	0.35 mm (28 Gauge)	Orange	TE/AD3735-008/01
Detachable	37 mm	0.45 mm (26 Gauge)	Green	TE/AD3745-008/01
Detachable	50 mm	0.45 mm (26 Gauge)	Blue	TE/AD5045-008/01
Detachable	75 mm	0.45 mm (26 Gauge)	Violet	TE/AD7545-008/01



**Packing:** 25 individually packed electrodes per box  
**Pre-wired:** 0.75 m

Reusable Cable	Wire Length	Part Number
Shielded	1.0 m	TE/D54-451
Unshielded	0.6 m	TE/D34-231



## Specifications

Not made with natural rubber latex	✓
REACH & RoHS (EU) compliant	✓
EtO sterilized	✓
Lead wire core	Tin Plated Copper
Lead wire insulation	PVC
Connector	DIN 42802-1
Shelf life for unopened pouches	36 months

Your Technomed Business Partner:



Legal manufacturer:

Technomed Europe

Wiebachstraat 25a

6466 NG Kerkrade

The Netherlands

T +31 43 408 68 68 | infotechnomed@welcony.com

Technomed USA

1626 Terrace Drive

Roseville | MN 55113 | United States

T +1 651 288 0890 | infotechnomed@welcony.com

www.technomedmedical.com