

Reusable EEG Cup Electrode Silver/Silver Chloride



Technomed understands the importance of a crystal clear signal, ease of use and electrode durability.

Technomed's Reusable Ag/AgCl EEG Cup Electrode is designed to provide a reliable signal and longtime ease of use.

The electrode head is designed with a larger dome volume for long term use.

The connector is optimized for easy headbox insertion and extraction.


User Benefits

- Clear and reliable signal with optimized layer of frequency stable Ag/AgCl on a solid Silver cup.
- Extra large dome volume for longer use compared to a standard dome.
- Short cup-to-cable connection to reduce lift risk.
- Optimal patient comfort with soft and flexible lead wires.
- Optimized connector for easy headbox insertion and extraction.
- Long electrode life with watertight design and Kevlar reinforced wires.
- Available in two different cup sizes.

Available variants

Materials	Cup Size	Cup Height	Wire Length	Part Number
Silver/Silver Chloride (Ag/AgCl)	10 mm	3.3 mm	1.0 m	TE/C22-434/99
Silver/Silver Chloride (Ag/AgCl)	10 mm	3.3 mm	1.5 m	TE/C22-634/99
Silver/Silver Chloride (Ag/AgCl)	10 mm	3.3 mm	2.0 m	TE/C22-834/99
Silver/Silver Chloride (Ag/AgCl)	10 mm	3.3 mm	2.5 m	TE/C22-934/99
Silver/Silver Chloride (Ag/AgCl)	6 mm	2.3 mm	1.5 m	TE/C24-634/99

Packing: sets of 12, 2 x 6 colors (red, green, white, yellow, blue, brown)



Specifications

Cup	Sterling Silver / AgCl layer
Not made with natural rubber latex	√
REACH & RoHS (EU) compliant	√
Lead wire design	Tinsel wire
Lead wire core	Kevlar and Silver Plated Copper
Lead wire insulation	TPU
Spring loaded connector	Gold Plated Brass / PVC
Connector	DIN 42802-1

Cleaning and disinfection: see IFU for details

Your Technomed Business Partner:



www.technomedmedical.com

Legal manufacturer::

Technomed Europe

Wiebachstraat 25a

6466 NG Kerkrade

The Netherlands

T +31 43 408 68 68 | infotechnomed@welcony.com

Technomed USA

1626 Terrace Drive

Roseville | MN 55113 | United States

T +1 651 288 0890 | infotechnomed@welcony.com